



## A Real-Time, Large-Area Transparent FLASH Beam Monitor and Dosimeter

Peter Friedman (PI & PD) – Integrated Sensors, LLC, Palm Beach Gardens, FL

Daniel Levin (PI), Nick Ristow, Claudio Ferretti, Monica Tecchio – University of Michigan Physics

Dale Litzenberg – University of Michigan Hospital, Ann Arbor, MI

Reinhard Schulte, Vladimir Bashkirov – Loma Linda University, Loma Linda, CA

P. S. Friedman, CEO | Integrated Sensors, LLC, Palm Beach Gardens, Florida, USA | www.isensors.net June 14, 2023



## Safety Problem: Monitoring FLASH Delivery

- To be clinically approved by the FDA/NIH/IRB, FLASH-RT must be monitored in real-time so the beam can be *immediately* terminated if a problem develops, given the <u>1000X\_higher dose rates</u> and ~ 0.1 to 0.5 sec total delivery time.
- Conventional beam monitors are not capable of large-area 2D imaging with sampling and analysis in real-time at FLASH rates (~ 50 µs).
- SOLUTION: Beam monitor that continuously images & analyzes FLASH-RT beam every 50 µs (20,000 fps) as the patient is irradiated.



# **FLASH-RT Monitor Performance Objectives**



#### $\checkmark$ = achieved

- Fast, real-time response to produce beam interlock, IEC compliant (in progress)
  - Deliver < 10% of total dose in case of fault
  - Scanning <u>proton</u> beams: process at 20 kfps → 50 µs, analysis in < 1 µs √</li>
  - Pulsed <u>electron</u> beams: data acquisition/process rate 100-1000 Hz → 1 ms (in progress)
- Spatial resolution: < 100  $\mu$ m on centroids  $\checkmark$
- Dynamic range: beam center to tails span 2 orders-of-magnitude ✓
- Low mass & thin profile: < 0.7 mm WE and profile/depth of ~ 11 cm  $\checkmark$
- Large area: 15 cm x 15 cm ✓ to 30 cm x 30 cm (in progress)
- Real-Time Dosimetry (within ~ 50  $\mu$ s for protons): < 4%  $\checkmark$
- Radiation resistance: 1 year of clinical usage  $(5d/wk) \le 1\%$  signal loss  $\checkmark$



## Large-Area FLASH Beam Scintillator Monitor

ullet



FBSM Conceptual Design

- Large 30 cm x 30 cm sensitive area
- Quadrant system with 4 cameras
- Thin 11 cm profile using folded optics
- Ultrafast machine-vision cameras
- Triggered or quasi-free-running modes
- FPGA data processing & analysis
- Low mass profile < 700 µm WE



## Large-Area Proprietary Scintillators

Type 1: Hybrid Material (HM) – Inorganic polycrystalline ceramic hybrid

• Thin < 700 µm WE

#### Type 2: Polymer Material (PM) – Semicrystalline

• Ultrathin to Thin: tested 2  $\mu$ m WE to < 300  $\mu$ m WE

#### Both Types 1 & 2 have favorable properties:

- Radiation hard
- Sharp images no internal reflections
- Non-hygroscopic
- Transmissive
- High light emittance for their respective type



## Test Beam Results from Prototypes

### FLASH Dose Rates

Notre Dame Radiation Lab ..... e<sup>-</sup> 8 MeV, 0.2-<u>235</u> Gy/s

U. Michigan Ion Beam Laboratory ..... p+ 5 MeV, > 300 Gy/s

Facility for Rare Isotope Beams ...... Kr<sup>+26</sup> 2.75 MeV/u, > 50 Gy/s

### Scintillator Efficiency Comparisons to Benchmarks



3 mm collimated electron beam ( $\beta^{-}$  source <sup>90</sup>Sr)

#### PM type





HM type



# Beam Current in HM Scintillator

#### Kr<sup>+26</sup> Beam Monitor Rate vs FRIB Instruments Rate



→ Result 1: Beam monitor measures currents over range covered by 4 different FRIB devices

- Faraday Cup
- Calibrated beam attenuator
- MCP detector
- Silicon detector

Result 2: Linear for 5 decades up to FLASH dose rates



## Dose Response in HM Scintillator

#### 8 MeV Electron Beam at 30 Hz (Notre Dame Lab)



Each data point = average of 100 pulses.

Pulse width in all cases = 2 ns.

At 30 Hz, average dose rate ranges from 6 to 234 Gy/s, which is well beyond 40 Gy/s for FLASH.

<u>Linear response</u> from Conventional-RT dose rate, to FLASH-RT high dose rate.



## Radiation Hardness of HM Scintillator



Low signal loss of 0.02%/kGy measured over 212 kGy.

> 1 yr of continuous clinical use≈ 1 % signal loss

Signal loss is correctable with internal UV calibration system.

# Beam Tests at U. Michigan Radiation Oncology





Varian linac: electron energy 6-16 MeV Conventional dose rate = 1-10 Gy/min



Protoype FLASH monitor single camera, 15 cm x 15 cm



# Beam Shape & Spatial Resolution

HM Scintillator vs Gafchromic Film



Beam monitor resolution is equal to Gafchromic film resolution (≲ 25 μm).

2D beam profiles are nearly identical.

Beam monitor primary advantage is real-time analysis

# 





We have demonstrated prototype monitors for FLASH-RT beams 🗸

- 2D Imaging with large area 15 cm x 15 cm  $\checkmark$  30 x 30 cm in development
- High sensitivity & dynamic range: single-particles to high FLASH level for RT
- Linear response for FLASH dose rates ✓
- Spatial resolution comparable with Gafchromic film: < 25 um  $\checkmark$
- Excellent radiation hardness
  - PM scintillator radiation damage: <u>none observed to 9 kGy</u> 🗸 Ο
  - HM scintillator radiation damage: <u>overall -0.02 %/kGy</u>, tested to 212 kGy 🗸 0
- Real-time FLASH data processing
  - $\circ$  20 kHz for protons with < 1  $\mu$ s required for data analysis  $\checkmark$
  - $\circ$  1 kHz electrons (in progress)

P. S. Friedman, CEO (peter@isensors.net) | Integrated Sensors, LLC, Palm Beach Gardens, Florida, USA | www.isensors.net







# FBSM Competitive Advantages

The FBSM is a patented, enabling technology for FLASH-RT

- <u>Two New High Efficiency Scintillators</u>
  - PM-scintillator (polymer) ultra-thin rolls
    HM-scintillator (hybrid) highest efficiency
- Innovative UFT Patented Configurations
  - <u>Ultra-fast</u> beam analysis ~ <u>50 µs</u> (continuous)
  - Real-time dosimetry, beam position & shape
  - > High spatial resolution (~ 25  $\mu$ m)
  - ➢ Water-equivalent thickness <700 µm</p>
  - Internal calibration
  - Multiple cameras & folded optics
  - > Detector area: up to  $\sim 30 \text{ cm x} 30 \text{ cm}$





# Applications

- **FLASH-RT** (electrons, protons, ions, x-rays, neutrons)
- Electron FLASH IORT (Intraoperative radiation therapy)
- Advanced EBRT including heavy-ions (alphas, carbon-ions, etc.)
- Beamline Monitoring (medical & scientific applications e.g., FRIB, ANL-ATLAS)
- Boron Neutron Capture Therapy (BNCT)
- Spatially Fractionated EBRT (Grid, Lattice, Minibeam, Microbeam)



# Validation & Collaborations

- Leading Academic & Government Institutions
  - University of Michigan
  - Loma Linda University
  - University of Texas / MD Anderson
  - Texas A&M
  - Notre Dame University
  - Florida State University
  - DOE Argonne National Laboratory
  - DOE Facility for Rare Isotope Beam





# Transmission in 6 µm PM Scintillator

<sup>86</sup>Kr<sup>+26</sup> ion (2.75 MeV/u) beam current measured with scintillator in/out of beam

